

Zinc finger protein GLI2

Catalog No: #AP83767



Package Size: #AP83767-1 50ug #AP83767-2 100ug #AP83767-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Zinc finger protein GLI2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Tax helper protein,Thp
Accession No.	Uniprot ID:Q0VGT2
Uniprot	Q0VGT2
GeneID	14633;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gli2

Note: This product is for in vitro research use only