

Glutamine synthetase

Catalog No: #AP83768



Package Size: #AP83768-1 50ug #AP83768-2 100ug #AP83768-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glutamine synthetase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	GS, Glutamate decarboxylase, Glutamate--ammonia ligase, GLNS
Accession No.	Uniprot ID:P46410
Uniprot	P46410
GeneID	396944;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GLUL

Note: This product is for in vitro research use only