

Hyaluronan-mediated motility receptor

Catalog No: #AP83769



Package Size: #AP83769-1 50ug #AP83769-2 100ug #AP83769-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hyaluronan-mediated motility receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Intracellular hyaluronic acid-binding protein, Receptor for hyaluronan-mediated motility, CD168, Ihabp, Rhamm
Accession No.	Uniprot ID:P97779
Uniprot	P97779
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Hmnr

Note: This product is for in vitro research use only