

Glucosidase 2 subunit beta

Catalog No: #AP83776



Package Size: #AP83776-1 50ug #AP83776-2 100ug #AP83776-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glucosidase 2 subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	80K-H protein,Glucosidase II subunit beta,Protein kinase C substrate 60.1 kDa protein heavy chain,PKC ζ ,Vacuolar system-associated protein 60,VASAP-60
Accession No.	Uniprot ID:Q28034
Uniprot	Q28034
GeneID	338067;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PRKCSH

Note: This product is for in vitro research use only