

Inhibin beta B chain

Catalog No: #AP83778



Package Size: #AP83778-1 50ug #AP83778-2 100ug #AP83778-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Inhibin beta B chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Activin beta-B chain
Accession No.	Uniprot ID:P17491
Uniprot	P17491
GeneID	25196;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Inhbb

Note: This product is for in vitro research use only