

Fatty acid desaturase 1

Catalog No: #AP83781



Package Size: #AP83781-1 50ug #AP83781-2 100ug #AP83781-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fatty acid desaturase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Delta(5) fatty acid desaturase,D5D
Accession No.	Uniprot ID:Q920L1
Uniprot	Q920L1
GeneID	76267;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fads1

Note: This product is for in vitro research use only