

Potassium/sodium hyperpolarization-activated cyclic nucleotide-gated channel 4

Catalog No: #AP83783

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP83783-1 50ug #AP83783-2 100ug #AP83783-3 1mg

Description

Product Name	Potassium/sodium hyperpolarization-activated cyclic nucleotide-gated channel 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Brain cyclic nucleotide-gated channel 3,BCNG-3,Bcng3
Accession No.	Uniprot ID:O70507
Uniprot	O70507
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Hcn4

Note: This product is for in vitro research use only