

78 kDa glucose-regulated protein

Catalog No: #AP83784



Package Size: #AP83784-1 50ug #AP83784-2 100ug #AP83784-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	78 kDa glucose-regulated protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	GRP-78,Heat shock 70 kDa protein 5,Immunoglobulin heavy chain-binding protein,BiP,GRP78
Accession No.	Uniprot ID:Q90593
Uniprot	Q90593
GeneID	396487;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HSPA5

Note: This product is for in vitro research use only