

Glutamate--cysteine ligase catalytic subunit

Catalog No: #AP83788



Package Size: #AP83788-1 50ug #AP83788-2 100ug #AP83788-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glutamate--cysteine ligase catalytic subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	GCS heavy chain, Gamma-ECS, Gamma-glutamylcysteine synthetase, GlcIc
Accession No.	Uniprot ID:P97494
Uniprot	P97494
GeneID	14629;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GlcIc

Note: This product is for in vitro research use only