

GSK3-beta interaction protein

Catalog No: #AP83790



Package Size: #AP83790-1 50ug #AP83790-2 100ug #AP83790-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	GSK3-beta interaction protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	GSKIP,RCJMB04_2i15
Accession No.	Uniprot ID:Q5ZMC6
Uniprot	Q5ZMC6
GeneID	423442;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GSKIP

Note: This product is for in vitro research use only