

Follistatin

Catalog No: #AP83793



Package Size: #AP83793-1 50ug #AP83793-2 100ug #AP83793-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Follistatin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FS,Activin-binding protein
Accession No.	Uniprot ID:P19883
Uniprot	P19883
GeneID	10468;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FST

Note: This product is for in vitro research use only