

Galactoside 3(4)-L-fucosyltransferase

Catalog No: #AP83795



Package Size: #AP83795-1 50ug #AP83795-2 100ug #AP83795-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Galactoside 3(4)-L-fucosyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Blood group Lewis alpha-4-fucosyltransferase, Lewis FT, FUTB, Fucosyltransferase 3, Fucosyltransferase III, FucT-III
Accession No.	Uniprot ID:Q11126
Uniprot	Q11126
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: FUT3

Note: This product is for in vitro research use only