

Insulin-like growth factor-binding protein 3

Catalog No: #AP83797



Package Size: #AP83797-1 50ug #AP83797-2 100ug #AP83797-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Insulin-like growth factor-binding protein 3 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | IBP-3,IBP3 |
| Accession No. | Uniprot ID:P17936 |
| Uniprot | P17936 |
| GeneID | 3486; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IGFBP3

Note: This product is for in vitro research use only