

Gasdermin-C

Catalog No: #AP83798



Package Size: #AP83798-1 50ug #AP83798-2 100ug #AP83798-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Gasdermin-C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Melanoma-derived leucine zipper-containing extranuclear factor, Mlze
Accession No.	Uniprot ID:P85967
Uniprot	P85967
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Gsdmc

Note: This product is for in vitro research use only