

# Procollagen galactosyltransferase 1

Catalog No: #AP83800



Package Size: #AP83800-1 50ug #AP83800-2 100ug #AP83800-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Procollagen galactosyltransferase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Collagen beta(1-O)galactosyltransferase 1,Glycosyltransferase 25 family member 1,Hydroxylysine galactosyltransferase 1,GLT25D1,PSEC0241
Accession No.	Uniprot ID:Q8NBJ5
Uniprot	Q8NBJ5
GeneID	79709;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:COLGALT1

Note: This product is for in vitro research use only