

# Solute carrier family 2, facilitated glucose transporter member 2



Catalog No: #AP83801

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP83801-1 50ug #AP83801-2 100ug #AP83801-3 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Solute carrier family 2, facilitated glucose transporter member 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Glucose transporter type 2, liver, GLUT-2, GLUT-2, GLUT2
Accession No.	Uniprot ID:Q90592
Uniprot	Q90592
GeneID	396272;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: SLC2A2

Note: This product is for in vitro research use only