

# Potassium/sodium hyperpolarization-activated cyclic nucleotide-gated channel 2

Catalog No: #AP83803

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP83803-1 50ug #AP83803-2 100ug #AP83803-3 1mg

## Description

Product Name	Potassium/sodium hyperpolarization-activated cyclic nucleotide-gated channel 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Brain cyclic nucleotide-gated channel 2,BCNG-2,BCNG2
Accession No.	Uniprot ID:Q9UL51
Uniprot	Q9UL51
GeneID	610;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:HCN2

Note: This product is for in vitro research use only