

Heparin cofactor 2

Catalog No: #AP83804



Package Size: #AP83804-1 50ug #AP83804-2 100ug #AP83804-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Heparin cofactor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Heparin cofactor II,HC-II,Protease inhibitor leuserpin-2,Serpin D1,Hcf2
Accession No.	Uniprot ID:Q64268
Uniprot	Q64268
GeneID	79224;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Serpind1

Note: This product is for in vitro research use only