

Glutaredoxin-1

Catalog No: #AP83816



Package Size: #AP83816-1 50ug #AP83816-2 100ug #AP83816-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Glutaredoxin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Thioltransferase-1,TTase-1,GRX
Accession No.	Uniprot ID:P12309
Uniprot	P12309
GeneID	397463;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GLRX

Note: This product is for in vitro research use only