

Hypoxanthine-guanine phosphoribosyltransferase

Catalog No: #AP83825



Package Size: #AP83825-1 50ug #AP83825-2 100ug #AP83825-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hypoxanthine-guanine phosphoribosyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	HGPRT,Hprt
Accession No.	Uniprot ID:P27605
Uniprot	P27605
GeneID	103689983;24465;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hprt1

Note: This product is for in vitro research use only