

Beta-hexosaminidase subunit beta

Catalog No: #AP83826

Package Size: #AP83826-1 50ug #AP83826-2 100ug #AP83826-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-hexosaminidase subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	65 kDa epididymal boar protein, Beta-N-acetylhexosaminidase subunit beta, Hexosaminidase subunit B, N-acetyl-beta-glucosaminidase subunit beta
Accession No.	Uniprot ID:Q29548
Uniprot	Q29548
GeneID	396958;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0) gene_full_name: HEXB

Note: This product is for in vitro research use only