

Hypoxia-inducible factor 1-alpha

Catalog No: #AP83827



Package Size: #AP83827-1 50ug #AP83827-2 100ug #AP83827-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hypoxia-inducible factor 1-alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	HIF-1-alpha
Accession No.	Uniprot ID:O35800
Uniprot	O35800
GeneID	29560;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Hif1a

Note: This product is for in vitro research use only