

Hepatocyte nuclear factor 1-beta

Catalog No: #AP83828



Package Size: #AP83828-1 50ug #AP83828-2 100ug #AP83828-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hepatocyte nuclear factor 1-beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	HNF-1-beta, Homeoprotein LFB3, Transcription factor 2, TCF-2, Variant hepatic nuclear factor 1, vHNF1, TCF2
Accession No.	Uniprot ID:P35680
Uniprot	P35680
GeneID	6928;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: HNF1B

Note: This product is for in vitro research use only