

Isopentenyl-diphosphate Delta-isomerase 1

Catalog No: #AP83830



Package Size: #AP83830-1 50ug #AP83830-2 100ug #AP83830-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Isopentenyl-diphosphate Delta-isomerase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Isopentenyl pyrophosphate isomerase 1,IPP isomerase 1
Accession No.	Uniprot ID:Q1LZ95
Uniprot	Q1LZ95
GeneID	514293;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IDI1

Note: This product is for in vitro research use only