

# Interferon-induced protein with tetratricopeptide repeats 1

Catalog No: #AP83831

Package Size: #AP83831-1 50ug #AP83831-2 100ug #AP83831-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interferon-induced protein with tetratricopeptide repeats 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IFIT-1, Interferon-induced 56 kDa protein, IFI-56K, G10P1, IFI56, IFNA1, ISG56
Accession No.	Uniprot ID:P09914
Uniprot	P09914
GeneID	3434;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: IFIT1

Note: This product is for in vitro research use only