

Immunoglobulin heavy constant gamma 1

Catalog No: #AP83832



Package Size: #AP83832-1 50ug #AP83832-2 100ug #AP83832-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Immunoglobulin heavy constant gamma 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Ig gamma-1 chain C region,Ig gamma-1 chain C region EU,Ig gamma-1 chain C region KOL,Ig gamma-1 chain C region NIE
Accession No.	Uniprot ID:P01857
Uniprot	P01857
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IGHG1

Note: This product is for in vitro research use only