

Gremlin-1

Catalog No: #AP83838



Package Size: #AP83838-1 50ug #AP83838-2 100ug #AP83838-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Gremlin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cell proliferation-inducing gene 2 protein,Cysteine knot superfamily 1, BMP antagonist 1,DAN domain family member 2,Down-regulated in Mos-transformed cells protein,Increased in high glucose protein 2,IHG-2,CKTSF1B1,DAND2,DRM,PIG2
Accession No.	Uniprot ID:O60565
Uniprot	O60565
GeneID	26585;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GREM1

Note: This product is for in vitro research use only