

Insulin-like growth factor-binding protein 6

Catalog No: #AP83839



Package Size: #AP83839-1 50ug #AP83839-2 100ug #AP83839-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Insulin-like growth factor-binding protein 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IBP-6,IBP6
Accession No.	Uniprot ID:P24592
Uniprot	P24592
GeneID	3489;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IGFBP6

Note: This product is for in vitro research use only