

# 3-hydroxyacyl-CoA dehydrogenase type-2

Catalog No: #AP83841



Package Size: #AP83841-1 50ug #AP83841-2 100ug #AP83841-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	3-hydroxyacyl-CoA dehydrogenase type-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	17-beta-hydroxysteroid dehydrogenase 10,17-beta-HSD 10,3-hydroxy-2-methylbutyryl-CoA dehydrogenase,3-hydroxyacyl-CoA dehydrogenase type II,Endoplasmic reticulum-associated amyloid beta-peptide-binding protein,Mitochondrial ribonuclease P protein 2,Mitochondrial RNase P protein 2,Type II HADH,Erab,Hadh2
Accession No.	Uniprot ID:O08756
Uniprot	O08756
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Hsd17b10

Note: This product is for in vitro research use only