

N-formyl peptide receptor 2

Catalog No: #AP83845



Package Size: #AP83845-1 50ug #AP83845-2 100ug #AP83845-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	N-formyl peptide receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FMLP-related receptor I,FMLP-R-I,Formyl peptide receptor-like 1,HM63,Lipoxin A4 receptor,LXA4 receptor,RFP,FPRH1,FPRL1,LXA4R
Accession No.	Uniprot ID:P25090
Uniprot	P25090
GeneID	2358;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FPR2

Note: This product is for in vitro research use only