

G-protein coupled estrogen receptor 1

Catalog No: #AP83849



Package Size: #AP83849-1 50ug #AP83849-2 100ug #AP83849-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	G-protein coupled estrogen receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Chemoattractant receptor-like 2,G protein-coupled estrogen receptor 1,G-protein coupled receptor 30,G-protein coupled receptor 41,Membrane estrogen receptor,mER,Cmkr12,Gper,Gpr30,Gpr41
Accession No.	Uniprot ID:O08878
Uniprot	O08878
GeneID	171104;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gper

Note: This product is for in vitro research use only