

HLA class I histocompatibility antigen, alpha chain E

Catalog No: #AP83854



Package Size: #AP83854-1 50ug #AP83854-2 100ug #AP83854-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HLA class I histocompatibility antigen, alpha chain E
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MHC class I antigen E,HLA-6.2,HLAE
Accession No.	Uniprot ID:P13747
Uniprot	P13747
GeneID	3133;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HLA-E

Note: This product is for in vitro research use only