

Heterogeneous nuclear ribonucleoprotein F

Catalog No: #AP83855



Package Size: #AP83855-1 50ug #AP83855-2 100ug #AP83855-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heterogeneous nuclear ribonucleoprotein F
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	hnRNP F, Nucleolin-like protein mcs94-1, HNRPF
Accession No.	Uniprot ID:P52597
Uniprot	P52597
GeneID	3185;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: HNRNPF

Note: This product is for in vitro research use only