

## Transcription factor HES-1

Catalog No: #AP83858



Package Size: #AP83858-1 50ug #AP83858-2 100ug #AP83858-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Transcription factor HES-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Class B basic helix-loop-helix protein 39,bHLHb39,Hairy and enhancer of split 1,Hairy homolog,Hairy-like protein,hHL,BHLHB39,HL,HRV
Accession No.	Uniprot ID:Q14469
Uniprot	Q14469
GeneID	3280;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:HES1

Note: This product is for in vitro research use only