

Hepcidin

Catalog No: #AP83861



Package Size: #AP83861-1 50ug #AP83861-2 100ug #AP83861-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Hepcidin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Hepc
Accession No.	Uniprot ID:Q99MH3
Uniprot	Q99MH3
GeneID	84604;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Hamp

Note: This product is for in vitro research use only