

Hepatocyte nuclear factor 1-alpha

Catalog No: #AP83863



Package Size: #AP83863-1 50ug #AP83863-2 100ug #AP83863-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Hepatocyte nuclear factor 1-alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	HNF-1-alpha,Transcription factor 1,TCF-1,TCF1
Accession No.	Uniprot ID:Q90867
Uniprot	Q90867
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HNF1A

Note: This product is for in vitro research use only