

Hematopoietic prostaglandin D synthase

Catalog No: #AP83865



Package Size: #AP83865-1 50ug #AP83865-2 100ug #AP83865-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hematopoietic prostaglandin D synthase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	H-PGDS,GST class-sigma,Glutathione S-transferase,Glutathione-dependent PGD synthase,Glutathione-requiring prostaglandin D synthase,Prostaglandin-H2 D-isomerase,Gsts,Pgds,Ptgds2
Accession No.	Uniprot ID:O35543
Uniprot	O35543
GeneID	58962;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hpgds

Note: This product is for in vitro research use only