

Bifunctional arginine demethylase and lysyl-hydroxylase JMJD6



Catalog No: #AP83867

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP83867-1 50ug #AP83867-2 100ug #AP83867-3 1mg

Description

Product Name	Bifunctional arginine demethylase and lysyl-hydroxylase JMJD6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Histone arginine demethylase JMJD6, JmjC domain-containing protein 6, Jumonji domain-containing protein 6, Lysyl-hydroxylase JMJD6, Peptide-lysine 5-dioxygenase JMJD6, Phosphatidylserine receptor, Protein PTDSR, PTDSR, RCJMB04_1m8
Accession No.	Uniprot ID:Q5ZMK5
Uniprot	Q5ZMK5
GeneID	417355;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: JMJD6

Note: This product is for in vitro research use only