

## Hepatoma-derived growth factor-related protein 3

Catalog No: #AP83870



Package Size: #AP83870-1 50ug #AP83870-2 100ug #AP83870-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Hepatoma-derived growth factor-related protein 3       |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Ms   |
| Other Names        | HRP-3  |
| Accession No.      | Uniprot ID:Q9JMG7                                      |
| Uniprot            | Q9JMG7   |
| GeneID             | 29877;   |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Hdgfrp3

Note: This product is for in vitro research use only