

3-hydroxy-3-methylglutaryl-coenzyme A reductase

Catalog No: #AP83871



Package Size: #AP83871-1 50ug #AP83871-2 100ug #AP83871-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	3-hydroxy-3-methylglutaryl-coenzyme A reductase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	HMG-CoA reductase
Accession No.	Uniprot ID:A7Z064
Uniprot	A7Z064
GeneID	407159;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: HMGCR

Note: This product is for in vitro research use only