

Glycerol kinase

Catalog No: #AP83878



Package Size: #AP83878-1 50ug #AP83878-2 100ug #AP83878-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glycerol kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	GK,ATP-stimulated glucocorticoid-receptor translocation promoter,ASTP,ATP:glycerol 3-phosphotransferase,Gyk
Accession No.	Uniprot ID:Q63060
Uniprot	Q63060
GeneID	79223;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gk

Note: This product is for in vitro research use only