

Interleukin-20 receptor subunit alpha

Catalog No: #AP83883



Package Size: #AP83883-1 50ug #AP83883-2 100ug #AP83883-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-20 receptor subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IL-20 receptor subunit alpha,IL-20R1
Accession No.	Uniprot ID:Q6PHB0
Uniprot	Q6PHB0
GeneID	237313;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL20ra

Note: This product is for in vitro research use only