

Integrin alpha-M

Catalog No: #AP83885



Package Size: #AP83885-1 50ug #AP83885-2 100ug #AP83885-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Integrin alpha-M
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CD11 antigen-like family member B,CR-3 alpha chain,Cell surface glycoprotein MAC-1 subunit alpha,Leukocyte adhesion receptor MO1,Neutrophil adherence receptor,CD11b,CD11B,CR3A
Accession No.	Uniprot ID:P11215
Uniprot	P11215
GeneID	3684;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ITGAM

Note: This product is for in vitro research use only