

Transcription factor 4

Catalog No: #AP83886

Package Size: #AP83886-1 50ug #AP83886-2 100ug #AP83886-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transcription factor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	TCF-4, Class A helix-loop-helix transcription factor ME2, Immunoglobulin transcription factor 2, ITF-2, SL3-3 enhancer factor 2, SEF-2, Itf2, Sef2
Accession No.	Uniprot ID:Q60722
Uniprot	Q60722
GeneID	21413;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0) gene_full_name: Tcf4

Note: This product is for in vitro research use only