

Golgi apparatus protein 1

Catalog No: #AP83887



Package Size: #AP83887-1 50ug #AP83887-2 100ug #AP83887-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Golgi apparatus protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CFR-1,Cysteine-rich fibroblast growth factor receptor,E-selectin ligand 1,ESL-1,Golgi sialoglycoprotein MG-160,CFR1,ESL1,MG160
Accession No.	Uniprot ID:Q92896
Uniprot	Q92896
GeneID	2734;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GLG1

Note: This product is for in vitro research use only