

Solute carrier family 2, facilitated glucose transporter member 3

Catalog No: #AP83892

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP83892-1 50ug #AP83892-2 100ug #AP83892-3 1mg

Description

Product Name	Solute carrier family 2, facilitated glucose transporter member 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Glucose transporter type 3, brain, GLUT-3, GLUT3
Accession No.	Uniprot ID:P11169
Uniprot	P11169
GeneID	6515;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SLC2A3

Note: This product is for in vitro research use only