

HLA class I histocompatibility antigen, alpha chain F

Catalog No: #AP83893



Package Size: #AP83893-1 50ug #AP83893-2 100ug #AP83893-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|----------------------------------------------------------------------------|
| Product Name | HLA class I histocompatibility antigen, alpha chain F |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | CDA12,HLA F antigen,Leukocyte antigen F,MHC class I antigen F,HLA-5.4,HLAF |
| Accession No. | Uniprot ID:P30511 |
| Uniprot | P30511 |
| GeneID | 3134; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HLA-F

Note: This product is for in vitro research use only