

## High mobility group protein HMG-I/HMG-Y

Catalog No: #AP83894



Package Size: #AP83894-1 50ug #AP83894-2 100ug #AP83894-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	High mobility group protein HMG-I/HMG-Y
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	HMG-I(Y),High mobility group AT-hook protein 1,High mobility group protein A1
Accession No.	Uniprot ID:Q8K585
Uniprot	Q8K585
GeneID	117062;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Hmga1

Note: This product is for in vitro research use only