

High mobility group protein B2

Catalog No: #AP83895



Package Size: #AP83895-1 50ug #AP83895-2 100ug #AP83895-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	High mobility group protein B2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	High mobility group protein 2,HMG-2,HMG2
Accession No.	Uniprot ID:P40673
Uniprot	P40673
GeneID	540444;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HMGB2

Note: This product is for in vitro research use only