

Growth hormone receptor

Catalog No: #AP83896



Package Size: #AP83896-1 50ug #AP83896-2 100ug #AP83896-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Growth hormone receptor |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Dog |
| Other Names | GH receptor,Somatotropin receptor |
| Accession No. | Uniprot ID:Q9TU69 |
| Uniprot | Q9TU69 |
| GeneID | 403721; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GHR

Note: This product is for in vitro research use only